

www.unitronusa.com

UNITRON[®] SINCE 1952

microscopes • components • video • binoculars • spotting scopes • telescopes • classic brass

Anywherescope™



This is a combination plastic and rubber armored microscope/telescope which can be used anywhere as the name suggests. It weighs approximately 23 ounces. The Anywherescope can truly be used anywhere and for a wide audience such as children and adults while camping, hiking, viewing coins, stamps, etc.



Features:

- 20x Microscope
- 7x Telescope
- Built in illuminator (uses a 3.0 volt lithium photo battery)
- Compass is built into the base
- The stage plate can be reversed and one side includes a cup to view liquids .
- The bottom of the base includes a ¼-20 thread that accepts any standard tripod or window mount.
- Vinyl case with carrying strap.
- The battery and bulb (2.47 volts lamp type #14 bulb type G-3 ½ average life is 15 hours) are easily changed by removing one thumb screw.

Much more information is available on our website (www.unitronusa.com) or by contacting Unitron directly.

This document is available on our website in Adobe® Acrobat® format as "[Unitron_Anywherescope_01.pdf](#)"

UNITRON[®]
S I N C E 1 9 5 2

Tel: (631) 589-6666 • Fax: (631) 589-6975

Email: info@unitronusa.com

All specifications subject to change without prior notice.
Prices are FOB Bohemia, New York, USA.

Unitron, Inc • 170 Wilbur Place • PO Box 469 • Bohemia, NY 11716