

www.unitronusa.com

UNITRON[®] SINCE 1952

microscopes • components • video • binoculars • spotting scopes • telescopes • classic brass




Biological Microscopes...

Unitron's complete line of specialty biological microscopes offers a model for every application — from very economical student level instrumentation to high-end laboratory units.

Our general line of biological microscopes include the following types, any of which may be augmented with our wide range of components...

- **MS Laboratory** — For biological and medical applications. Fulfills medical students' requirements at realistic prices. Available in monocular, binocular and trinocular models.



MS Laboratory

- **Polarizing Laboratory** — Performance and versatility usually found on more expensive models. Includes strain-free objectives, built-in bertrand lens, compensators, graduated rotatable stage and micrometer eyepiece.



Biozoom

- **Biozoom** — Biological zoom microscopes with magnifications with one objective from 50X to 300X. Available in either laboratory or field models.



Polarizing Laboratory

- **Unilab Laboratory** (not shown) — Routine medical or college classroom laboratory microscope. Choice of achromatic or planachromatic objectives.

Much more information is available on our website (www.unitronusa.com) or by contacting Unitron directly.

This document is available on our website in Adobe® Acrobat® format as "[Unitron_Biologicals_01.pdf](#)"

UNITRON[®]
S I N C E 1 9 5 2

Tel: (631) 589-6666 • Fax: (631) 589-6975

Email: info@unitronusa.com

All specifications subject to change without prior notice.
Prices are FOB Bohemia, New York, USA.

Unitron, Inc • 170 Wilbur Place • PO Box 469 • Bohemia, NY 11716