


www.unitronusa.com

UNITRON[®] SINCE 1952



microscopes • components • video • binoculars • spotting scopes • telescopes • classic brass

Metallographs...

For nearly 50 years, the recognized reputation for dependable high quality metallographs and metallurgical microscopes has made UNITRON the specified choice in leading industrial, educational and research facilities throughout the world. Our products are easily adaptable to video, camera or computer.



- **VERSAMET** - Research metallograph with unique optical system providing instant sequential selection of four illumination techniques brightfield, darkfield, polarized light and Nomarski interference contrast. Built-in camera system with fully automatic photo exposure and CCTV adapter.
- **NEOMET** - Inverted metallurgical microscope offers brightfield and polarized light. Designed for ease of operation for both visual photographic and CCTV usage.
- **MS2-IM2** - Inverted metallurgical microscope for basic applications. Accessories for photography and CCTV.
- **MS2 and MS2 BD** (not shown) - Upright materials microscope for routine applications. Monocular, binocular and trinocular. Two models: MS2 Brightfield and MS2 BD Brightfield, Darkfield and Polarized Light.
- **RMM** - Portable depth measuring rollscope measuring microscope with X,Y and Z motion. Ideal for observing/photographing surfaces too large for ordinary microscopes.



Much more information is available on our website (www.unitronusa.com) or by contacting Unitron directly.

This document is available on our website in Adobe® Acrobat® format as "[Unitron_Metallographs_01.pdf](#)"

UNITRON[®]
S I N C E 1 9 5 2

Tel: (631) 589-6666 • Fax: (631) 589-6975

Email: info@unitronusa.com

All specifications subject to change without prior notice.
Prices are FOB Bohemia, New York, USA.

Unitron, Inc • 170 Wilbur Place • PO Box 469 • Bohemia, NY 11716