

www.unitronusa.com

UNITRON[®] SINCE 1952

microscopes • components • video • binoculars • spotting scopes • telescopes • classic brass



Instruments for Special Applications...

FORENSIC PROJECTORS - MP6B Measuring and CP-6 Comparison (not shown) for measuring and comparison of forensic material such as bullets, cartridges, fibers, paint, etc.

SHOPSCOPE - Practical, all purpose modular shop inspection microscope. Wide variety of interchangeable components to meet inspection and measuring needs. Magnification ranges 10x to 200x.



Surface Measuring



MP6B Projector

MASK ALIGNMENT MICROSCOPES - Specifically designed for manufacturers of instruments used for the accurate alignment of masks or parts in the microcircuit industry. These microscopes give a split field of view offering models with 24mm to 100mm, 5mm to 30mm and 5mm to 80mm separation between objectives.



Shopscope



Hole Inspection Microscope

HOLE INSPECTION MICROSCOPE - Special optics for inspecting wall surfaces of through hole diameters .002" to .0375". Up to 6" deep. Zoom feature gives full field of view. Typical applications include circuit boards, spinnerettes and internal plating.

SURFACE MEASURING INSTRUMENTATION - Compact surface measuring utilizing the stylus method that includes unique touch screen multi-layer control screen that displays functions and parameters necessary for accurate measurements.



Mask Alignment Microscopes



Zoom Objective

Much more information is available on our website (www.unitronusa.com) or by contacting Unitron directly.

This document is available on our website in Adobe® Acrobat® format as "[Unitron_SpecialApps_01.pdf](#)"

UNITRON[®]

SINCE 1952

Tel: (631) 589-6666 • Fax: (631) 589-6975

Email: info@unitronusa.com

All specifications subject to change without prior notice.
Prices are FOB Bohemia, New York, USA.

Unitron, Inc • 170 Wilbur Place • PO Box 469 • Bohemia, NY 11716