

www.unitronusa.com

UNITRON[®] SINCE 1952

microscopes • components • video • binoculars • spotting scopes • telescopes • classic brass

Stereo Microscopes and Accessories...

Long a standard in the industry, Unitron's line of stereo microscopes provide a wide selection to meet customer needs.

These stereo microscopes include the following, any of which may be augmented with our wide range of components...

- **ZF SERIES STEREO** - Convergent prism optics for best resolution, quality and field flatness includes fixed and zoom models with magnification ranges 7x to 180x. Trinocular and binocular available.
- **M SERIES STEREO** - Parallel beam system offers the greatest depth of focus. Includes one zoom model, three fixed models and three dual discussion models - magnification ranges from 6x to 200x.



Triad of Stereo Microscopes
Continuous variable zoom magnification



The ZF Series



The M Series



Remember - Unitron supplies a complete line of microscope accessories, too!

Much more information is available on our website (www.unitronusa.com) or by contacting Unitron directly.

This document is available on our website in Adobe[®] Acrobat[®] format as "[Unitron_Stereo_01.pdf](#)"

All specifications subject to change without prior notice.
Prices are FOB Bohemia, New York, USA.

UNITRON[®]
S I N C E 1 9 5 2

Tel: (631) 589-6666 • Fax: (631) 589-6975

Email: info@unitronusa.com

Unitron, Inc • 170 Wilbur Place • PO Box 469 • Bohemia, NY 11716