

www.unitronusa.com

UNITRON[®] SINCE 1952

microscopes • components • video • binoculars • spotting scopes • telescopes • classic brass

MS Stereo Microscopes

MSZSB

The Unitron MSZSB 0.7x - 4.5x Zoom Stereo Microscope features a parfocal twin zooming prismatic objective system with the capability of magnifications of 3.5x to 180x with auxiliary lenses and eyepieces. Assembly and inspection in industry with the flexibility of greater magnification and the ability of rotating the head. Ideal for hospitals, medical facilities, teaching, stamp/coin collecting, wire bonding, agriculture, forestry, geology and archeology.

MSZSM

The Unitron MSZSM 0.75x - 3.4x Zoom Stereo Microscope features a parfocal twin zooming prismatic objective system. Applications include: Assembling and inspection in industry. Teaching in schools to view materials in the sciences. Viewing gross materials in hospitals and medical facilities. Storing images in a computer for future reference. Other uses include agriculture, forestry, geology, archaeology, and museums.

MSFSB

The Unitron MSFSB Fixed Binocular Stereo Microscopes are available in 6 models. The MSFSB are high quality prism systems with parfocal objectives which give versatility where there is only an occasional need to change objectives. Modular system offers flexibility. Ideal for assembly and inspection where only one or two magnifications are needed. Least expensive type of stereo microscope make them popular for schools, field work, and to check water samples for contamination. Examining rocks, bugs, and other gross samples.



UNITRON[®]
MS

UNITRON[®]

SINCE 1952

Tel: (631) 589-6666 • Fax: (631) 589-6975

Email: info@unitronusa.com

Much more information is available on our website (www.unitronusa.com) or by contacting Unitron directly.

This document is available on our website in Adobe® Acrobat® format as "Unitron_Stereo_02.pdf"

All specifications subject to change without prior notice.
Prices are FOB Bohemia, New York, USA.

Unitron, Inc • 170 Wilbur Place • PO Box 469 • Bohemia, NY 11716

