

PRESS RELEASE

FOR IMMEDIATE RELEASE

Contact: *Mark Clymer*
Director of Marketing
UNITRON Ltd.
267-733-6725

UNITRON Introduces the Z12 Zoom Stereo Microscope

Commack, NY [February 15, 2018] — — UNITRON Ltd. is pleased to introduce the Z12 zoom stereo microscope. With a magnification range from 8X to 400X and three configurable and adjustable bases, the Z12 addresses a wide variety of applications from assembly and quality control in the industrial markets, to stem cell and embryo selection in the life sciences. The new UNITRON Z12 delivers both performance and high-end features at a value-conscious price point.

The Z12 includes many features typically found in higher priced microscopes. First, the Z12 is available with three new stands to suit the application: plain focusing stand, transmitted light base with sliding/tilting mirror, and a slim profile LED transmitted light base. Color-balanced LED light sources give bright illumination to the specimen, and are highly cost effective (low energy consumption and long lifespans). A 5°-45° tilting trinocular observation head offers ergonomic viewing comfort that is customizable to the observer, and apochromatic objectives deliver crisp color reproduction and high resolution of fine structures. Click-stop zoom positions allow the observer to quickly return to a specific magnification setting, thus saving time and increasing reproducibility. “Our customers have asked for these features that are usually found only in premium branded microscopes,” says Brian Taub, Executive Vice President of UNITRON. “We are delighted to give them those features, high performance and more in the new UNITRON Z12.”

For more information on the UNITRON Z12 zoom stereo microscope, please visit
<http://www.unitronusa.com/products/stereo-microscopes/z12-zoom-stereo-microscope-series>.

####

About UNITRON

UNITRON specializes in high quality, precision instruments for industrial, clinical, life science, research and educational applications. Since 1952 UNITRON's microscopes and related optical accessories have been consistently used and trusted worldwide in companies such as Intel, BAE Systems, ITT, GE, Williams Co., DuPont, Boeing, 3M, Alcoa, MIT, Raytheon, Sony, Texas Instruments, iRobot, Northrop Grumman, United Technologies, and the Mayo Clinic. Manufactured and assembled to the strictest quality standards, our products are available exclusively through a nationwide network of authorized distributors. Visit us at <http://unitronusa.com/>.