

Z730 ZOOM STEREO MICROSCOPE SERIES

QUICK REFERENCE GUIDE

The UNITRON Z730 Series with a 4.3:1 zoom ratio is an ideal choice for industrial inspection, assembly, quality control, OEM integration, research and education applications. Featuring outstanding optics and a large field of view with a high depth of focus, the Z730 reduces sample handling times allowing for faster inspections and a more efficient work process.

Viewing Heads



45° Inclined
Binocular Head
#132-11-45



60° Inclined
Binocular Head
#132-11-60



45° Inclined
Trinocular Head
#132-11-70

Viewing Heads

Item#	Description
132-11-45	Inclined 45° Binocular
132-11-60	Inclined 60° Binocular
132-11-70	Inclined 45° Trinocular

Eyepieces



WF10x/20mm FOV
#132-10-10



WF15x/16mm FOV
#131-10-15



WF20x/12.5mm FOV
#131-10-20

Eyepieces

Item#	Description
132-10-10	10x/20mm Eyepiece
131-10-15	15x/16mm Eyepiece
131-10-20	20x/12.5mm Eyepiece

Objectives



0.5x/W D 180mm
Auxiliary
#132-15-05



0.63x/W D 140mm
Auxiliary
#132-15-06



0.75x/W D 120mm
Auxiliary
#132-15-75



2.0x / W D 30mm
Auxiliary
#132-15-02

Auxiliary Objectives

Item#	Description
132-15-05	0.5x (Use with #131-20-40)
132-15-06	0.63x
132-15-75	0.75x
132-15-02	2.0x

Adapters



0.5x C-mount
#132-25-05



0.5x Auxiliary Lens
Riser Adapter
#131-20-40

Adapters

Item#	Description
132-25-05	0.5x C-mount Adapter
131-20-40	Riser Adapter (Use with #132-15-05)