




15854 LED RING LIGHT

- 72 LED bulbs with 10 step variable intensity light control
- Sectional lighting control; allows users to direct the light from 4 separate quadrants
- 60.5mm inner diameter
- Color Temperature: 6400 Kelvin
- Lamp life: 10,000 hours
- Power supply: AC 90—265v, 50/60Hz
- Output: 12v 400mA
- RoHS compliant, ESD safe



15852 LED RING LIGHT

- 60 LED bulbs with variable light control and external power supply.
- 62mm inner diameter
- Lamp life: 10,000 hours
- Power supply: AC 90—265v, 50/60Hz
- Output: 24v 6w
- External power supply dimensions: 122 x 94 x 44 mm
- ESD safe




15853 POLARIZING LED RING LIGHT

- Variable LED ring light with removable polarizer attachment.
- The polarizer reduces glare from highly reflective surfaces and is easily rotatable.
- 70 mm inner diameter
- Lamp life: 10,000 hours
- Power supply: AC 90—265v, 50/60Hz
- Output: 24v 6w
- External power supply dimensions: 122 x 94 x 44 mm
- ESD safe



15856 LED RING LIGHT

- 60 LED bulbs with variable intensity light control
- ESD safe ABS plastic housing
- 60mm inner diameter
- Lamp life: 10,000 hours
- Power supply: AC 90 - 265v, 50/60Hz
- Output: 6 watt
- ESD safe



15857 LED DUAL GOOSENECK ILLUMINATOR

- Dual gooseneck variable intensity LED illuminator
- Each gooseneck has a 3 watt LED bulb allowing for pinpoint light control
- Lamp life: 10,000 hours
- Power supply: AC 100-240v, 50/60Hz
- ESD safe

UNITRON®

73 Mall Drive, Commack, New York 11725

Tel: (631) 543-2000 Fax: (631) 589-6975

Email: info@unitronusa.com • Website: www.unitronusa.com

**ISO 9001
Certification**

Design and production
adheres to ISO9001
international quality standard.

**ISO 14001
Certification**

Design and production meets the
requirements of international standard
ISO 14001 for environmental management.