

# AU-500-MADM - MINI AUTOMATIC DIGITAL MICROSCOPE



A good assistant for engineers  
An amazing tool for the curious

- Industrial Inspection
- Quality Control
- Dermatology
- Education
- Art Preservation
- Hobbies & Collecting




# Features — AU-500-MADM

Magnification	12x - 200x
Frame Rates	VGA - Up to 40fps 720P - Up to 30fps 5.0MP - Up to 15fps
Focus Mode	Auto Focus
Color Balance	Auto White Balance Auto Exposure Auto Gain
Illumination	Built-in LEDs
Interface	High-speed USB 2.0
Dimensions	Φ50.5 mm x 45 mm
Weight	75g

*Design, features and specifications are subject to change without notice.*

## Simple Push Button Controls

- Zoom in/out
- Autofocus
- Brightness level
- Snapshot



### Auto Focus

The built-in auto focus lens and the push button feature makes image focusing easy!

### Built-in LEDs

Adjust the bright white LED lights to 16%, 33%, 50%, 66%, 83% or 100%!



### Recommended PC Specs:



- 1.2 GHz or higher
- 512MB RAM
- 20GB or higher hard drive
- CD-ROM drive for installation
- Windows XP, 2000, Vista, 7, Windows Server 2003

### Built-in Microphone

Allows you to record an analysis notation!



Auto-focus



White lights



5.0MP



Hi-Speed USB 2.0



Take pictures



Microphone



Remote synchronization



Measurement and analysis



Video recording



Real color restore



View and preview



12x~200x continuously zoom



#### ISO 9001 Certification

Design and production adheres to ISO9001 international quality standard.

#### ISO 14001 Certification

Design and production meets the requirements of international standard ISO 14001 for environmental management.



73 Mall Drive • Commack, NY 11725  
631-864-1000 (P) • 631-543-8900 (F)  
info@accu-scope.com • info@unitronusa.com  
www.accu-scope.com • www.unitronusa.com

