

Z650HR ZOOM STEREO MICROSCOPE SERIES

QUICK REFERENCE GUIDE

With a built-in camera port and a 0.55x c-mount adapter (0/100 light split), cameras can be directly mounted to the Unitron Z650HR. Featuring an 8:1 zoom ratio (0.6x – 5.0x), a 115mm working distance and WF10x/23mm focusable eyepieces, the Z650HR is ideal for industrial inspection applications.

Viewing Head



30° Inclined
Trinocular Head
#134-11-03

Viewing Heads

Item#	Description
134-11-03	Inclined Trinocular With built-in 0.55x c-mount adapter

Eyepieces



WF10x/23mm FOV
#134-10-11



WF15x/16mm FOV
#131-10-15



WF20x/12.5mm FOV
#131-10-20



WF30x/8mm FOV
#134-10-30

Eyepieces

Item#	Description
134-10-11	10x/23mm Eyepiece
131-10-15	15x/16mm Eyepiece
131-10-20	20x/12.5mm Eyepiece
134-10-30	30x/8mm Eyepiece

Objectives



0.5x/WD 211mm
Auxiliary
#131-15-05



0.7x/WD 137mm
Auxiliary
#131-15-07



2.0x/WD 43.5mm
Auxiliary
#131-15-02

Auxiliary Objectives

Item#	Description
131-15-05	0.5x
131-15-07	0.7x
131-15-02	2.0x

Adapters & Accessories



Fiber Optic
Ring Light Adapter
#131-15-00



Polarizer/Analyzer
#131-30-00

Adapters & Accessories

Item#	Description
131-15-00	Fiber Optic Ring Light Adapter
131-30-00	Polarizer/Analyzer